

# **HSAPath: Activation Analysis for Human Cellular Pathways**

## Addendum to ESKAPE Act PLUS Manual

This application extends the functionality of ESKAPE Act PLUS to the human genome. Given data from high-throughput experiments from human subjects, HSAPath returns information on activated and repressed KEGG pathways and GO terms.

Most of manual developed for ESKAPE Act PLUS applies to HSAPath. The sections that are unique to the ESKAPE Act PLUS version of the application are highlighted in yellow. Section V, “Worked Example”, outlines an analysis using differential gene expression data from pseudomonas samples. There is updated sample data for HSAPath, which can be downloaded in the application via the “HSA example data” button. After downloading the data you can follow the same procedure outlined in Section V, ignoring the steps regarding selection of bacterial strain.

The data file containing KEGG pathways, GO terms and their constituent genes was constructed using the EnrichmentBrowser (version 2.28.0) and the limma (version 3.54.1) R packages. Enrichment Browser was used to get GO terms and limma was used to get KEGG pathways.



# ESKAPE Act **PLUS**

**Activation Analysis for ESKAPE Pathogens  
and other **P**rokaryotes **L**abs **U**sually **S**tudy**

User Manual

Second edition: September 2022

Thomas Hampton, PhD  
Department of Microbiology and Immunology  
The Geisel School of Medicine at Dartmouth  
Hanover, NH 03755, USA  
e-mail: [Thomas.H.Hampton@Dartmouth.edu](mailto:Thomas.H.Hampton@Dartmouth.edu)

**If using this app advances your research, please cite our publication:**

[ESKAPE Act Plus: Pathway Activation Analysis for Bacterial Pathogens.](#)

Koeppen K, Hampton TH, Neff SL, Stanton BA. *mSystems*. 2022 Oct 19:e0046822.

doi: [10.1128/msystems.00468-22](https://doi.org/10.1128/msystems.00468-22). PMID: 36259735.

## Contents

<b>I. Quick Start Guide .....</b>	<b>3</b>
1. Input file upload .....	3
2. <b>Strain selection</b> .....	3
3. Run activation analysis.....	3
4. Obtain results .....	3
5. Reset application.....	3
<b>II. Introduction.....</b>	<b>4</b>
1. ESKAPE pathogens .....	4
2. KEGG pathways and GO terms.....	4
3. Activation analysis versus enrichment analysis .....	4
4. Why use ESKAPE Act PLUS? .....	5
<b>III. Input requirements .....</b>	<b>6</b>
1. Input file format .....	6
2. UniProt ID conversion .....	6
<b>IV. Overview of capabilities .....</b>	<b>10</b>
1. <b>Supported strains</b> .....	10
2. Number of KEGG pathways and GO terms for each strain .....	10
3. Statistics .....	11
4. App output.....	12
<b>V. Worked example .....</b>	<b>13</b>
<b>VI. Troubleshooting .....</b>	<b>18</b>
<b>VII. References.....</b>	<b>19</b>

## I. Quick Start Guide

### 1. Input file upload

Click the “Browse” button to upload a csv file with input data. The file should contain UniProt stable entry identifiers in column 1 and log2 fold changes between two treatment groups in column 2.

Example data can be downloaded by clicking the “Pseudomonas PA14 example data” button.

### 2. Strain selection

From the dropdown menu, select the strain that corresponds to the UniProt entry IDs in the input file.

When using the example data, choose “Pseudomonas aeruginosa PA14”.

### 3. Run activation analysis

Once an input file has been uploaded, a “Run Activation Analysis” button will appear. Click on it to start the Analysis.

This user manual can be downloaded with the “User Manual” button.

### 4. Obtain results

Results will appear in tabs in the main panel of the application (more details on page 12). Output files can be downloaded with the “Download results” button.

### 5. Reset application


The application can be reset to its original state with the “Reset” button.

**Upload CSV file:**

pau\_Result.csv

Upload complete

Required file format: csv file needs to contain UniProt stable entry identifier in column 1 and log2 fold changes between two treatment groups in column 2.

 Pseudomonas PA14 example data


File has the correct number of columns.


**Select strain:**

Pseudomonas aeruginosa PA14 ▼

Upload a csv file with UniProt protein identifiers and log2 fold changes, select a strain, then click on 'Run Activation Analysis' that will appear below.

Run Activation Analysis

 User Manual

 Download results

Reset

## II. Introduction

### 1. ESKAPE pathogens

The group of bacteria collectively known as ESKAPE (*Enterococcus faecium*, *Staphylococcus aureus*, *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Pseudomonas aeruginosa* and *Enterobacter sp.*) are pathogens of concern due to their virulence and ability to develop antibiotic resistance. Collectively, they are the major cause of life-threatening hospital-acquired infections world-wide. The ESKAPE Act PLUS web application enables analysis of high-throughput experiments on ESKAPE pathogens and several other strains of bacteria that are of relevance to the biomedical research community. A complete list of supported strains is provided on page 10.

### 2. KEGG pathways and GO terms

ESKAPE Act PLUS uses biological pathway information from the [Kyoto Encyclopedia of Genes and Genomes](#) (KEGG) (M. Kanehisa and Goto 2000), obtained through the KEGGREST R package version 1.32.0 (Tenenbaum 2021), and [gene ontology](#) (GO) term annotations (Ashburner et al. 2000; Gene Ontology Consortium 2021) for biological processes, molecular functions and cellular components retrieved from the [Universal Protein Resource](#) (UniProt) (The UniProt Consortium 2021). KEGG pathway visualizations in the app are made possible via an API interface with KEGG Mapper's [Color Tool](#) (Minoru Kanehisa and Sato 2020; Minoru Kanehisa, Sato, and Kawashima 2022). In these visualizations, KEGG pathway images are overlaid with gene or protein level fold changes, allowing for easy visualization of pathway-level activation or repression. Genes or proteins with a positive fold change are highlighted in magenta, while genes or proteins with a negative fold change are depicted in green. Genes or proteins that were not detected in the experiment are shown in white. ESKAPE Act PLUS provides links to downloadable pathway level images for all significant KEGG pathways.

### 3. Activation analysis versus enrichment analysis

Many publicly available tools allow users to perform a gene set enrichment analysis for a biological pathway or GO term. Enrichment analysis typically requires a criterion such as a threshold p-value to identify a subset of differentially expressed genes or proteins. Differentially expressed genes are said to be “enriched” in genes performing specific biological functions if the subset of differentially expressed genes contains a higher proportion of genes performing these functions than one would expect by chance. Enrichment therefore suggests an association between experimental conditions and biological function. However, selecting different p-value cutoffs can

change enrichment results, and typical choices such as  $FDR < 0.05$  may result in very small gene sets, limiting power. In addition, enrichment analysis does not consider the direction of the fold changes and thus does not predict the biological effect associated with the induction or repression of a given pathway or GO term. By contrast, ESKAPE Act PLUS uses all genes or proteins that were detected in a high-throughput experiment and their associated fold changes to predict overall activation or repression at the level of a biological pathway or GO term based on the fold changes of all genes associated with a given pathway or GO term.

#### 4. Why use ESKAPE Act PLUS?

ESKAPE Act PLUS makes it easier for researchers working with prokaryotes to interpret high-throughput experiments such as RNA-seq and proteomics, leading to new insights and hypotheses. Specifically, ESKAPE Act PLUS provides biological pathway and GO term activation analysis, which were previously only available for eukaryotic systems and popular model organisms.

### III. Input requirements

#### 1. Input file format

Input files need to be files with a comma separated values format (.csv) that contain two columns without headers:


- **Column 1** contains UniProt identifiers (ID conversion is explained in section 2. below).
- **Column 2** contains log2 fold changes for the UniProt identifier specified in Column 1, for the comparison of interest. **This column must contain numeric values only, no text, and may not contain any missing values.**

Selecting “Download Pseudomonas PA14 example data” will download an example input data file for *Pseudomonas aeruginosa* strain PA14 (“pau\_Result.csv”) to illustrate the correct data format.

#### 2. UniProt ID conversion

Step 1 – go to the [UniProt homepage](#)

Step 2 – select “ID mapping”:



The screenshot shows the UniProt homepage. The navigation bar at the top includes links for BLAST, Align, Peptide search, ID mapping (highlighted with a red dashed circle), and SPARQL. The main header reads "Find your protein" with a search bar and a dropdown menu set to "UniProtKB". Below the search bar, there are examples: "Insulin, APP, Human, P05067, organism\_id:9606". A footer section contains four tiles: "Proteins UniProt Knowledgebase" (Reviewed: 568,002, Unreviewed: 226,771,948), "Species Proteomes" (Protein sets for species with sequenced genomes from across the tree of life), "Protein Clusters UniRef" (Clusters of protein sequences at 100%, 90% & 50% identity), and "Sequence Archive UniParc" (Non-redundant archive of publicly available protein sequences seen across different databases).

**Step 3** – paste your gene identifiers into the provided window or upload a file containing your gene identifiers:

## Retrieve/ID mapping

Enter your IDs or [load from a text file](#). Separate IDs by whitespace (space, tab, newline) or commas.

ABAYE2456  
 ABAYE2463  
 ABAYE2436  
 ABAYE2418  
 ABAYE2460  
 ABAYE2415  
 ABAYE2444  
 ABAYE2416  
 ABAYE2465  
 ABAYE2419  
 ABAYE2450  
 ABAYE2438  
 ABAYE2442  
 ABAYE2445

ⓘ Your input contains 3750 IDs

**Step 4** – Select from “Gene name” (or other identifier, if applicable) to “UniProtKB”, enter the Taxon ID for your strain of interest, select the organism once the popup window appears and click the “Map IDs” button. A table with the Taxon IDs for all supported strains can be found on page 10. The example below is for *Acinetobacter baumannii* stain AYE (Taxon ID 509173):

From database  
 Gene Name

To database  
 UniProtKB

Restrict by taxonomy  
 509173  
 Acinetobacter baumannii (strain AYE) [509173]

Name your ID Mapping job  
 ABAYE2456 +3749 Gene\_Name → UniProtKB

Reset Map 3750 IDs

**Step 5** – On the following page, select the most recent ID mapping job:

## Tool results

Your tool analysis results from the last ⓘ 7 days are listed below. If you have tools jobs running, you can navigate away to other pages and you will be notified once the job is completed.

Job type	Name	Created	Status	
ID MAPPING	ABAYE2456 +3749 Gene_Na	2022-08-25 10:19	Completed (607 hits)	☆ 📄 🗑
c101d6c44cdf0dd5aa92f1896f29177f491d790a				



**Step 6** – Verify that a sufficient number of identifiers were successfully mapped. Small numbers of mapped identifiers (e.g., fewer than 200) may reduce the number of significant KEGG pathways and GO terms that can be identified. Bear this in mind, then click the “Download” button:

Taxonomy

Filter by taxonomy

Status

Reviewed (Swiss-Prot) (310)

Unreviewed (TrEMBL) (3,237)

Proteins with

3D structure (15)

Active site (256)

Activity regulation (7)

Beta strand (3)

Binding site (425)

More items

Protein existence

Predicted (2,069)

Homology (1,461)

Protein level (16)

Transcript level (1)

Annotation score

4 (18)

3 (291)

2 (663)

1 (2,575)

ID mapping 3,607 results found for Gene\_Name → UniProtKB

Overview

Input Parameters

API Request

BLAST

Align

Map

Download

Add

View: Cards

Table

Customize columns

Share

143 IDs were not mapped:

Show IDs

From	Entry	Entry Name	Protein Names	Gene Names	Organism	Length
<input type="checkbox"/>	ABAYE2456	B0VAX6	B0VAX6_ACIBY	Beta-lactamase[...]	ABAYE2456	Acinetobacter baumannii (strain AYE) 424 AA
<input type="checkbox"/>	ABAYE2463	B0VAV9	B0VAV9_ACIBY	Uncharacterized protein	ABAYE2463	Acinetobacter baumannii (strain AYE) 543 AA
<input type="checkbox"/>	ABAYE2436	B0VAY0	B0VAY0_ACIBY	Uncharacterized protein	ABAYE2436	Acinetobacter baumannii (strain AYE) 302 AA
<input type="checkbox"/>	ABAYE2418	B0VB07	B0VB07_ACIBY	Uncharacterized protein	ABAYE2418	Acinetobacter baumannii (strain AYE) 133 AA
<input type="checkbox"/>	ABAYE2460	B0VAV7	B0VAV7_ACIBY	Putative hydroxyacyl-CoA dehydrogenase[...]	ABAYE2460	Acinetobacter baumannii (strain AYE) 316 AA
<input type="checkbox"/>	ABAYE2415	B0VB05	B0VB05_ACIBY	Uncharacterized protein	ABAYE2415	Acinetobacter baumannii (strain AYE) 167 AA
<input type="checkbox"/>	ABAYE2444	B0VAY8	B0VAY8_ACIBY	Uncharacterized protein	ABAYE2444	Acinetobacter baumannii (strain AYE) 105 AA
<input type="checkbox"/>	ABAYE2416	B0VB06	B0VB06_ACIBY	Uncharacterized protein	ABAYE2416	Acinetobacter baumannii (strain AYE) 177 AA

Feedback

Help

**Step 7** – select “TSV” from the dropdown menu and click “Download”:

**Download** [x]

☐ Download selected (0)  
☒ Download all (3,607)

Format  
TSV

Compressed  
☒ Yes  
☐ No

Customize columns

Reviewed x Entry Name x Protein names x Gene Names x Organism x Length x

Data 6 External links

Search for available columns [q]

Names & Taxonomy	4	▼
Sequences	1	▼
Function		▼
Miscellaneous	1	▼
Interaction		▼
Expression		▼
Gene Ontology (GO)		▼
Pathology & Biotech		▼
Subcellular location		▼
PTM / Processing		▼

Generate URL for API Preview 10 Cancel **Download**

**Step 8** – Create the input data file for ESKAPE Act PLUS by merging the Tab-separated file downloaded from UniProt (containing the original gene identifiers and the UniProt entry IDs required by ESKAPE Act PLUS) with the table containing gene identifiers and fold changes. Instructions on how to merge data using different platforms can be found in the links below:

[Merging tables in Excel](#)

[Merging tables using R](#)

[Merging tables using Python](#)

## IV. Overview of capabilities

### 1. Supported strains

ESKAPE Act PLUS currently supports 23 strains from 13 species, including the six ESKAPE pathogens as well as several other species that are of relevance to the biomedical research community. Listed below are the KEGG, UniProt and Taxon identifiers for each strain.

Species	Strain	KEGG ID	UniProt ID	Taxon ID
<i>Acinetobacter baumannii</i>	AYE	aby	ACIBY	509173
<i>Acinetobacter baumannii</i>	MDR-ZJ06	abz	ACIBM	497978
<i>Bacteroides fragilis</i>	NCTC 9343	bfs	BACFN	272559
<i>Bacteroides ovatus</i>	3725 D1 iv	boa	BACOV	28116
<i>Bacteroides thetaiotaomicron</i>	VPI-5482	bth	BACTN	226186
<i>Clostridioides difficile</i>	630	cdf	CLOD6	272563
<i>Clostridioides difficile</i>	R20291	cdl	CLODR	645463
<i>Enterobacter cloacae</i>	ATCC 13047	enc	ENTCC	716541
<i>Enterococcus faecium</i>	DO	efu	ENTFD	333849
<i>Escherichia coli</i>	K-12 MG1655	eco	ECOLI	511145
<i>Escherichia coli</i>	K-12 W3110	ecj	ECOLI	316407
<i>Escherichia coli</i>	O157:H7 EDL933	ece	ECO57	155864
<i>Escherichia coli</i>	O157:H7 Sakai	ecs	ECO57	386585
<i>Escherichia coli</i>	BL21(DE3)	ebd	ECOBD	469008
<i>Klebsiella pneumoniae</i>	MGH 78578	kpn	KLEP7	272620
<i>Prevotella melaninogenica</i>	ATCC 25845	pmz	PREMB	553174
<i>Pseudomonas aeruginosa</i>	UCBPP-PA14	pau	PSEAB	208963
<i>Pseudomonas aeruginosa</i>	PAO1	pae	PSEAE	208964
<i>Staphylococcus aureus</i>	COL	sac	STAAC	93062
<i>Staphylococcus aureus</i>	USA300 FPR3757	saa	STAA3	367830
<i>Staphylococcus aureus</i>	NCTC8325	sao	STAA8	93061
<i>Staphylococcus aureus</i>	Newman	sae	STAAE	426430
<i>Streptococcus sanguinis</i>	SK36	ssa	STRSV	388919

### 2. Number of KEGG pathways and GO terms for each strain

ESKAPE Act PLUS will calculate activation for KEGG pathways and GO terms that have 4 or more associated genes in the input data provided. For example, if the input data contain only three genes on a specific path, no activation score will be calculated.

The total number of possible KEGG pathways and GO terms available for each strain are listed below.

Species	Strain	KEGG paths	GO terms
<i>Acinetobacter baumannii</i>	AYE	92	1622
<i>Acinetobacter baumannii</i>	MDR-ZJ06	93	1516
<i>Bacteroides fragilis</i>	NCTC 9343	85	1398
<i>Bacteroides ovatus</i>	3725 D1 iv	80	1285
<i>Bacteroides thetaiotaomicron</i>	VPI-5482	83	1482
<i>Clostridioides difficile</i>	630	83	1339
<i>Clostridioides difficile</i>	R20291	84	1240
<i>Enterobacter cloacae</i>	ATCC 13047	94	2210
<i>Enterococcus faecium</i>	DO	93	2310
<i>Escherichia coli</i>	K-12 MG1655	94	3936
<i>Escherichia coli</i>	K-12 W3110	94	3939
<i>Escherichia coli</i>	O157:H7 EDL933	94	2267
<i>Escherichia coli</i>	O157:H7 Sakai	70	1237
<i>Escherichia coli</i>	BL21(DE3)	92	2105
<i>Klebsiella pneumoniae</i>	MGH 78578	99	2164
<i>Prevotella melaninogenica</i>	ATCC 25845	99	2929
<i>Pseudomonas aeruginosa</i>	UCBPP-PA14	101	1921
<i>Pseudomonas aeruginosa</i>	PAO1	71	985
<i>Staphylococcus aureus</i>	COL	88	1378
<i>Staphylococcus aureus</i>	USA300 FPR3757	86	1357
<i>Staphylococcus aureus</i>	NCTC8325	85	1334
<i>Staphylococcus aureus</i>	Newman	82	1776
<i>Streptococcus sanguinis</i>	SK36	70	1191

### 3. Statistics

ESKAPE Act PLUS uses a binomial test to assess whether a given KEGG pathway or GO term contains significantly more induced or repressed genes or proteins than would be expected by chance. KEGG pathway and GO term activation analysis were performed as previously described (Hampton et al. 2018; Goodale et al. 2019; Koeppen et al. 2021). Under the null hypothesis, about 50% of the genes or proteins associated with any KEGG pathway or GO term would be expected to respond to treatment with a FC > 0, while the other 50% would respond with a FC < 0. Statistically significant divergence from this 50:50 percent split is assessed using binomial tests based on the fold changes of all genes or proteins associated with a KEGG pathway or GO term. The universe for the binomial test is all genes or proteins from the input file (rather than all genes or proteins known for a given species). Using identifiers found in the source file as the reference avoids the possibility of introducing biases related to tissue type and sample preparation.

Pathway enrichment p-values from the binomial test are FDR-corrected to account for the total number of KEGG pathways or GO terms that were tested. Only KEGG pathways or GO terms with an FDR-corrected p-value  $< 0.05$  are considered significantly induced or repressed and displayed in figures and tables of the in-app results output. The downloadable results .csv files also contain the uncorrected p-values for each KEGG pathway and GO term to allow users to identify borderline cases with uncorrected p-values  $< 0.05$  that just barely missed the significance cutoff of  $\text{FDR} < 0.05$ .

#### 4. App output

ESKAPE Act PLUS generates two kinds of output – in-app output and downloadable output. Following an activation analysis, the in-app output appears in tabs in the main panel of the application, while downloadable output becomes available via the “Download results” button that appears in the left control panel of the app.

##### **In-app output:**

- Number of user-provided UniProt Entry identifiers that mapped to the internal KEGG and GO database as well as the number of significant KEGG pathways and GO terms.
- Boxplots for all significant ( $\text{FDR} < 0.05$ ) KEGG pathways and GO terms.
- Tables with significant ( $\text{FDR} < 0.05$ ) KEGG pathways and GO terms, containing estimates ( $< 0.5$  = repressed,  $> 0.5$  = activated), median fold changes, p-values from the binomial test and FDR-corrected p-values.
- Pathway level visualizations: the KEGG table contains clickable links to KEGG Mapper pathway level visualizations for each significant KEGG pathway. See page 4 for more details on pathway level visualizations.

##### **Downloadable output:**

- Boxplots for all significant ( $\text{FDR} < 0.05$ ) KEGG pathways and GO terms.
- Tables with estimates ( $< 0.5$  = repressed,  $> 0.5$  = activated), median fold changes, p-values, and FDR-corrected p-values for all KEGG pathways and GO terms that were tested.
- html file (“Pathways.html”) with links to pathway level visualizations for all significant KEGG pathways generated by KEGG Mapper.

## V. Worked example

 Pseudomonas PA14 example data

1. Download *Pseudomonas aeruginosa* strain PA14 example data file “pau\_Result.csv”

### Upload CSV file:

Browse... pau\_Result.csv

Upload complete

Required file format: csv file needs to contain UniProt stable entry identifier in column 1 and log2 fold changes between two treatment groups in column 2.

 Pseudomonas PA14 example data

File has the correct number of columns.

### Select strain:

Pseudomonas aeruginosa PA14 ▼

Upload a csv file with UniProt protein identifiers and log2 fold changes, select a strain, then click on 'Run Activation Analysis' that will appear below.


Run Activation Analysis

2. Upload example data file “pau\_Result.csv” using the “Browse...” button.


3. Select strain “*Pseudomonas aeruginosa* PA14”

4. Click “Run Activation Analysis”

5. The **“Summary”** tab shows how many user-provided identifiers mapped to the internal database for KEGG pathway and GO term activation analysis and how many KEGG pathways and GO terms were significantly induced or repressed:



6. The **“KEGG Graph”** tab contains a boxplot of log<sub>2</sub> fold changes for each gene or protein on significant KEGG pathways. Magenta indicates pathway induction and green indicates repression:




7. The **“KEGG Table”** tab provides a searchable and sortable table of significant KEGG pathways with pathway activation estimates ( $< 0.5$  = repressed,  $> 0.5$  = activated), median fold changes, p-values from the binomial test and FDR-corrected p-values:

Summary **KEGG Graph** **KEGG Table** GO Graph GO Table

Show **50** entries Search:


Path	Estimate	Median_fold_change	P_value	FDR
<a href="#">Biofilm formation</a>	0.33	-0.1	0.00219490890778059	0.00949158066743457
<a href="#">Pyruvate metabolism</a>	0.69	0.07	0.00379370563537163	0.0132077159157383
<a href="#">Glyoxylate and dicarboxylate metabolism</a>	0.7	0.06	0.00456153274640824	0.0147856578676681
<a href="#">Alanine, aspartate and glutamate metabolism</a>	0.71	0.03	0.0166738478001207	0.0460982850944515

The Path column contains clickable links to pathway images generated by KEGG mapper, such as this example image:





8. The **“GO Graph”** tab contains a boxplot of log2 fold changes for each gene or protein associated with a significant GO term. Magenta indicates induction and green indicates repression:



9. The **“GO Table”** tab provides a searchable and sortable table of significant GO terms with pathway activation estimates ( $< 0.5$  = repressed,  $> 0.5$  = activated), median fold changes, p-values from the binomial test and FDR-corrected p-values:

Summary   KEGG Graph   KEGG Table   GO Graph   **GO Table**

Show  entries   Search:

GO_term	Estimate	Median_fold_change	P_value	FDR
DNA binding	0.61	0.02	0.000407810912196372	0.00649948641312968
integral component of membrane	0.62	0.04	0.00000341509554582281	0.000124407052026402
DNA-binding transcription factor activity	0.64	0.03	0.000142346183659903	0.00302485640277294


**Upload CSV file:**

Browse...

pau\_Result.csv

Upload complete

Required file format: csv file needs to contain UniProt stable entry identifier in column 1 and log2 fold changes between two treatment groups in column 2.


 Pseudomonas PA14 example data


File has the correct number of columns.


**Select strain:**

Pseudomonas aeruginosa PA14 ▼

Upload a csv file with UniProt protein identifiers and log2 fold changes, select a strain, then click on 'Run Activation Analysis' that will appear below.

 Download results

 Reset

 User Manual

10. All results generated by the app (boxplots of significant KEGG pathways and GO terms as well as csv files with results for all KEGG pathways and GO terms and an html file with links to KEGG Mapper pathway level visualizations) can be downloaded with the **“Download results”** button.

11. The **“Reset”** button returns the app to its initial state before running a new analysis.

## VI. Troubleshooting

Error message	Probable cause	Solution
Error – file has fewer than 2 columns!	The input file that was uploaded contains fewer than 2 columns of data	Upload an input file that contains UniProt entry IDs in column 1 and log2 fold changes in column 2
Error – file has more than 2 columns!	The input file that was uploaded contains more than 2 columns of data	Upload an input file that only contains 2 columns - UniProt entry IDs in column 1 and log2 fold changes in column 2
<b>ID error</b> - double check selected strain	The user-provided identifiers do not match the UniProt identifiers of the selected strain. Either the selected strain does not match the UniProt IDs of the input file or the input file does not contain the correct UniProt Entry IDs.	Double-check that the correct strain is selected that corresponds to the UniProt Entry IDs provided in the input file. If the input file does not contain the correct UniProt IDs, please refer to the section on ID conversion (page 6) to fix this problem.
<b>Data error</b> - The second column of the uploaded file needs to contain numeric values only, no missing values.	Verify that the second column of the input file contains a numeric fold change for each observation and does not contain any missing values.	Upload a corrected input file that contains a numeric value for each observation in the second column (no empty fields or rows).
<b>Data error</b> - The second column of the uploaded file needs to contain numeric values only, no text.	The second column of the input file contains text (rather than numeric values) for at least one of the rows. The input file contains a header with column names.	Upload a corrected input file without column name headers that contains numeric values only in the second column.

## VII. References

- Ashburner, M., C. A. Ball, J. A. Blake, D. Botstein, H. Butler, J. M. Cherry, A. P. Davis, et al. 2000. "Gene Ontology: Tool for the Unification of Biology. The Gene Ontology Consortium." *Nature Genetics* 25 (1): 25–29. <https://doi.org/10.1038/75556>.
- Gene Ontology Consortium. 2021. "The Gene Ontology Resource: Enriching a Gold Mine." *Nucleic Acids Research* 49 (D1): D325–34. <https://doi.org/10.1093/nar/gkaa1113>.
- Goodale, Britton C., Thomas H. Hampton, Emily N. Ford, Craig E. Jackson, Joseph R. Shaw, Bruce A. Stanton, and Benjamin L. King. 2019. "Profiling MicroRNA Expression in Atlantic Killifish (*Fundulus Heteroclitus*) Gill and Responses to Arsenic and Hyperosmotic Stress." *Aquatic Toxicology (Amsterdam, Netherlands)* 206 (January): 142–53. <https://doi.org/10.1016/j.aquatox.2018.11.009>.
- Hampton, Thomas H., Craig Jackson, Dawoon Jung, Celia Y. Chen, Stephen P. Glaholt, Bruce A. Stanton, John K. Colbourne, and Joseph R. Shaw. 2018. "Arsenic Reduces Gene Expression Response to Changing Salinity in Killifish." *Environmental Science & Technology* 52 (15): 8811–21. <https://doi.org/10.1021/acs.est.8b01550>.
- Kanehisa, M., and S. Goto. 2000. "KEGG: Kyoto Encyclopedia of Genes and Genomes." *Nucleic Acids Research* 28 (1): 27–30. <https://doi.org/10.1093/nar/28.1.27>.
- Kanehisa, Minoru, and Yoko Sato. 2020. "KEGG Mapper for Inferring Cellular Functions from Protein Sequences." *Protein Science: A Publication of the Protein Society* 29 (1): 28–35. <https://doi.org/10.1002/pro.3711>.
- Kanehisa, Minoru, Yoko Sato, and Masayuki Kawashima. 2022. "KEGG Mapping Tools for Uncovering Hidden Features in Biological Data." *Protein Science: A Publication of the Protein Society* 31 (1): 47–53. <https://doi.org/10.1002/pro.4172>.
- Koeppen, Katja, Amanda Nymon, Roxanna Barnaby, Laura Bashor, Zhongyou Li, Thomas H. Hampton, Amanda E. Liefeld, et al. 2021. "Let-7b-5p in Vesicles Secreted by Human Airway Cells Reduces Biofilm Formation and Increases Antibiotic Sensitivity of *P. Aeruginosa*." *Proceedings of the National Academy of Sciences of the United States of America* 118 (28): e2105370118. <https://doi.org/10.1073/pnas.2105370118>.
- Tenenbaum, Dan. 2021. *KEGGREST: Client-Side REST Access to the Kyoto Encyclopedia of Genes and Genomes (KEGG)* (version 1.32.0).
- The UniProt Consortium. 2021. "UniProt: The Universal Protein Knowledgebase in 2021." *Nucleic Acids Research* 49 (D1): D480–89. <https://doi.org/10.1093/nar/gkaa1100>.